



Recombinant Arabidopsis thaliana BTB/POZ domain-containing protein At4g30940 (At4g30940)

Product Code	CSB-MP525379DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O65555
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGLSNDRIKF NVGGRIFETT ATTLANAGRD SFFGALFDEN WNLSQPGDLF IDRNPDCFAV LLDLLRTGDL NIPPNIPEL LHKEAMFYGL IDHLRTAKWG PFDGNRLHLS RSVTGIAPGD GTAIRAGPDG GCCIAHGSVV HVFDWML EEH PTINLDYQRV NDVGWIDSGN IVLSACERLG RGDGGMGLFS SSSGELRYKF QVSHDNQVKS YSAGALSFSP DSKIFTSCKG RSNEYGIGVW DQSTGKQVDF FYESPGWSLG DADKLQWLSG KNCLLVATLF PRKDNCYISL LDFREKNMVW SWSDIGFLTM AEEKRVRDAI AMEESNSICV VNEFEDLGFI DLRMDGGGSS VRWSSRSRLM KSKMPDEPCY PKLALHEGQL FSSMNDISISV FCGSDWVLTS RLKRSYGGSI CDFSIGGDRL FALHSEENVF DVWETLPPPI I
Source	Mammalian cell
Target Names	At4g30940
Protein Names	Recommended name: BTB/POZ domain-containing protein At4g30940
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.