



Recombinant *Schizosaccharomyces pombe* Uncharacterized aminotransferase C777.03c (SPCC777.03c)

Product Code	CSB-YP525590SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74542
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTVSEGIARI QAFRKQVPLV EKNENITFFN LSFQPPMNSI VANTINEYVD WGLYNPNPKP LWQQKAEETR ELVAKYLNAS SSDSIAFTRD TTEGLNLFQR SMHFQPGDNV VLLDGEHPNH GFGWISLQEA GLEVRLIPSK DIYYADASTF EPYVDEKTKA IGISSVMFHS GQLNNVKDIC NKFRPENIHV LVDMTQQVGL SSIDVQELNV SACAFACHKG LSCPTGLGVL YVAPDVL PQL KNPPIV GAG AIENLDANLL VNLNPIYFST ARFEHLNKS LITTLCLNNY LSFLMDIGFE NIEQYLRELG KSLVDELSKL GVTIIGSKDP KKRSTHSYVA KILNPEWDAF LKNEGVYASQ YRDGIRLSLG LYNNAADISR LANTIKKGV L NKIPFN
Source	Yeast
Target Names	SPCC777.03c
Protein Names	Recommended name: Uncharacterized aminotransferase C777.03c EC= 2.6.-.-
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.