



# Recombinant Schizosaccharomyces pombe Peptide-N (4)- (N-acetyl-beta-glucosaminyl)asparagine amidase

<b>Product Code</b>	CSB-BP525600SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O74739
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDFHAISQRF IDMMRSKNSQ NASQPPETYP FYHEVRQMSQ HPWMYEDPEL QDYALSILPL DKLFQDASEL EKEGDGSWGY QDYVIQALLK WFKREFFVWV NQPPCEKCGG ETHMTGNGPP NEEEKWNGVR NVELYQCNVC GHNQRFPRYN RIRALLDSRK GRCGEWANCF TFLCRALGSR ARWIWNAEDH VWTEVYSNKQ QRWVHVDSGE ESFDEPLIYE QGWGKKMSYC LGFGIDSVRD VSHRYIRHPE NGLPRDRCPE SVLQQALHEI NIEFRSRLTD SERKALEEED KREKDELDDGY MRPVSQATPT NTDLPARQTG NVEWKEKRGE AGK
<b>Source</b>	Baculovirus
<b>Target Names</b>	png1
<b>Protein Names</b>	Recommended name: Peptide-N(4)-(N-acetyl-beta-glucosaminyl)asparagine amidase Short name= PNGase EC= 3.5.1.52 Alternative name(s): Peptide:N-glycanase 1
<b>Expression Region</b>	1-333
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.