



Recombinant Schizosaccharomyces pombe Putative box C/D snoRNA protein SPCC613.07 (SPCC613.07)

Product Code	CSB-YP525624SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74906
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSENRLGICS TCQKNASKYR CPRCDSRFCC LECNLEHKRL TKCSGERDPA TFVPSKSLVN HLNSDFNFLS GVERLINRKE NGSHQVSNRA ERNRLQLKRS LERAGINIKF APPSNKKRRL NRTHYDKKSH LIKWSIEWCL HESSTSKDLT DEASENTIIT HSHPESEPLE KIFRKLVEEN SEMNSQVSKA NIHAPDDMQF MIKSRKSRSK GFIIYKKIEPS NSLSSCLRNS FVFEVPTIHV FTSTTQVHTE SSYETSSSEQ SDDSSSSSSC ISDSEESSSD DELNELSNEK KANSPTQNV D TSSKLSSVKF QNEDDEDRRK NSDGSEASYS LPPLFGSYFQ KQGQY
Source	Yeast
Target Names	SPCC613.07
Protein Names	Recommended name: Putative box C/D snoRNA protein SPCC613.07
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.