



Recombinant Schizosaccharomyces pombe Peptidase M20 domain-containing protein C757.05c (SPCC757.05c)

Product Code	CSB-BP525626SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74916
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	GPLNF GYLSWGQHSI LEQDSVQLIS HQSLLPAGDN LISLHKSLVE IESLGGNEVN VSDFLKSYLE SKGLTVELQR VSSNPTVRDN VYAYLGSQRN TKVVLTSID TVNPFLPYYI EGDKIHGGRGS CDAKSSVAAQ IFAMLELMSE GKVQEGDLSL LFVVGEEIDG IGMKTVANKL NADWEVAIFG EPTENKLGVG HKGNFRFDVY AHGKACHSGY PQEGFSAIEF LLRQCVKLM DTLPKSKLLG PSTINIGTIE GGAAANILAA EAKAEVFIRV AEDIETIRDI AEDLFDTEHS EIKVIQYSPP QYLDYDIPGM DTVVLAYATD IPYLSDRKMK IYLFPGGSIR EAHGPNEYVT FSQLEGLNG YKRMVMYNLR
Source	Baculovirus
Target Names	SPCC757.05c
Protein Names	Recommended name: Peptidase M20 domain-containing protein C757.05c
Expression Region	26-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.