



Recombinant Human Jerky protein homolog (JRK)

Product Code	CSB-BP525639HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O75564
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSTLYDIRAH KAQLLRFFAS SDSNKALEQR RTLHTPKLEH LDRVLYEWFL GKRSEGVVPS GPMLIEKAKD FYEQMQLTEP CVFSGGWLWR FKARHGKIKL DASSEKQSAD HQAAEQFCFAF FRSLAAEHGL SAEQVYNAD E TGLFWRCLPN PTPEGGAVPG PKQGKDRLTV LMCANATGSH RLKPLAIGKC SGPRAFKGIQ HLPVAYKAQG NAWVDKEIFS DWFHHIFVPS VREHFRTIGL PEDSKAVLLL DSSRAHPQEA ELVSSNVFTI FLPASVASLV QPMEQGIRRD FMRNFINPPV PLQGPHARYN MNDAIFSVAC AWNVAVPSHFV RRAWRKLWPS VAFAEGSSSE EELEAEFCFPV KPHNKSFahi LELVKEGSSC PGQLRQRQAA SWGvAGREAE GGRPPAATSP AEVvWSSEKT PKADQDGRGD PGEgEEVAWE QAAVAFDaVL RFAERQPCFS AQEVGQLRAL RAVFRSQQQV RRRRGALGAV VKVEALQEGP GGCGATAQSP LPCSSTAGDN
Source	Baculovirus
Target Names	JRK
Protein Names	Recommended name: Jerky protein homolog
Expression Region	1-520
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.