



Recombinant *Drosophila atripex* Alpha-amylase-related protein (Amyrel)

Product Code	CSB-MP525650DKZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O77011
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila atripex</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	QHNPHWWGNR NTIVHLFEWK WSDIAAECES FLGPRGFAGV QVSPVNENII SAGRPWWERY QPISYKLVTR SGNEQEFADM VRRCDNVGVR IYVDVLLNHM SGDFDGI AVG TAGSEAEPSK KSYPGVPYSA LDFHPTCEIT DWNDRFQVQQ CELVGLKDLQ QSSEWVRSKL IEFLDHLIEL GVAGFRVDAK KHMAADDLSY IYSSISDLNI EHGFPHNARP FIFQEVIDHG HETVSREEYN QLGAVTEFRF SEEIGNAFRG NNALKWLQSW GTGWGFLPSG QALTFVDNHD NQRDMGAVLN YKSPKQYKMA TAFHLAYPYG ISRVMSFAF DDHDTAPPQD EQERIISPEF DEEGACVNGW ICEHRWRQIY AMVGFKNAVR DTELSNWWDN GDNQISFCRG NKGFLAVNNN LYDLSRELQT CLPAGVYCDV ISGSLIDGSC TGKSVTV DGN GYGYIHIGSD DFDGVLALHV DARI
Source	Mammalian cell
Target Names	Amyrel
Protein Names	Recommended name: Alpha-amylase-related protein EC= 3.2.1.1
Expression Region	21-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.