



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At2g29860 (At2g29860)

Product Code	CSB-MP525761DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82370
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVLLSEIPGG SNGDDPNMNP QELPEELIES IVAPIPCYY PSLSLLSRAF RHVITSQQLF VTRSGLGFKV PVLAFYFICT PYTTPRWFIL RRSNIPLQLR RLNSLPHMFP GAAVVTIGYK IYVMGGYNYQ PVSTVIIIDC RFHTWHYLQD MQRARYHATP GVIDGRIYVI GGRKKQDADW VEVFDVTTES WETVPTQCPN EASENFLVFT YAVMQGRILG VFLLTNQDKV YGSHWELYVS
Source	Mammalian cell
Target Names	At2g29860
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At2g29860
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.