



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At2g29800 (At2g29800)

Product Code	CSB-EP525762DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82376
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MASSETSDD GSNGGVNPQK PEEPHKNPQE EKENQENPN EAEDEDDHQD EEVENVPIIP RKIPVLIEN TIAPLRRCHY PKLSLLSNAF RQVISSDLF QVRSLIGSTE PVLTYLITFK YPTFEEGRWF ILQRRNNTSL KLNCVTSLPP MFLGCTAVTI GHKIYVVGYY NFRYNTIST VLEIDCRFNT CRHLRNMKRD RCSAVAGVID GRIYVVAGRQ RRFDDWVEVF DVETERWELV PGPFSSFASS SGKFIVHVVL DNKIYIMDGD YCFAYDPRRR RWETWGPESA QRSYWHLSSC VVDDLLYAIV PREIFGASIV VYDPRGIWR PVMGLEFWPN LVYFESKMAN FGGKLVILGC YRSSFDYYRK DVWCVEVALE KHEDGQIWGK VESLSLVNAF PMSPPFELSR TVTI
Source	E.coli
Target Names	At2g29800
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At2g29800
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.