



Recombinant Salmonella typhimurium Malonyl CoA-acyl carrier protein transacylase (fabD)

Product Code	CSB-YP525958SXB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O85140
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	TQFAFVFPG QGSQSVGMLA EMAANYPIVE ETFAEASAAL GYDLWALTQQ GPAEELNKTW QTQPALLTAS VALWRVWQQQ GGKMPALMAG HSLGEYSALV CAGVINFADA VRLVEMRGKF MQEAVPEGTG GMSAIIGLDD ASIAKACEES AEGQVVSPVN FNSPGQVVIA GHKEAVERAG AACKAAGAKR ALPLPVSVPVPS HCALMKPAAD KLAVELAKIT FSAPTVPVNV NVDVKCETDA AAIRDALVRQ LYNPVQWTKS VEFIAAQGVE HLYEVGPGKV LTGLTKRIVD TLTASALNEP AALSAALTQ
Source	Yeast
Target Names	fabD
Protein Names	Recommended name: Malonyl CoA-acyl carrier protein transacylase Short name= MCT EC= 2.3.1.39
Expression Region	2-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.