



Recombinant Fowlpox virus Putative ankyrin repeat protein FPV031 (ANK3)

Product Code	CSB-BP526032FEY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O90760
Product Type	Recombinant Protein
Immunogen Species	Fowlpox virus (strain NVSL) (FPV)
Purity	>85% (SDS-PAGE)
Sequence	MPFLIELLRI TGIVTLCPKY RLDNRNSLLLV ATKRNYIDVV RYLVKKGVDI NFQETIRDNL TPLMIASRFN SHQLVELLLN NGAIINQRSL TCGNTALHLA VKNDNRITVD ILLFHGANTN ITNNDGFTPL HKAVIYNASI DIIKLLRYK ADVNIIRDNEE ENTGLTPLDI AMSCNNYEII SLLVSHVIRL DYSTCMSNKT KGFDPHNKLV KDNKRLQRIA IHCLKDIEKM KQVSINSRFT LFDLFVNNNV DLLLRCINKD NRFIVNFDKK LTVFNILYTD FIDTFMSRHV LLNKASKVLE DLFLDNDSNT SSWNNLPNEI KDHFITYINN DELKIMTGSK T
Source	Baculovirus
Target Names	ANK3
Protein Names	Recommended name: Putative ankyrin repeat protein FPV031
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.