



Recombinant Schizosaccharomyces pombe Putative G3BP-like protein (nxt3)

Product Code	CSB-YP526082SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94260
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTAENATLLE PVLGKDEIGW MFVQEYYTYL NKEPNRLHCF YTKKSTLIHG DEGESISLCH GQQEIH NKIL DLDFQNCKVL ISNVDSLASS NGGIVIQVLG EMSNKGKLSR KFAQTFFLAE QPNGYFVLND IFRFLREDVE EEEESPDAVE KEKKDVASEP YVNGVQSQEH LPSAKEEGHY QDPAATENNF ATAALISNET DSL NQATLAV PEEPVIQVTE ASVPSFVSQQ ENQLQDEALT SNSKNADAIG ASDANVATAP KSWADLIARN HPDVKSQASV SSTASTTGQT VKGVNADQTQ QPTAPYTQSN ELLETSVFKV NIPPETSDVS LKSAMSIFGP VKAIEFARRK GTAYVDFVNH ECVQLALNKK TLQINNATLN IEERRRLFSG KFNKSGDKKS NDNYNGMKRN FRKGNRGAFD GRSKEVTTTSK KQNN
Source	Yeast
Target Names	nxt3
Protein Names	Recommended name: Putative G3BP-like protein
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.