



Recombinant Schizosaccharomyces pombe U3 small nucleolar RNA-associated protein 15 (utp15)

Product Code	CSB-MP526102SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94365
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSAPAGRLQP QRFAALPNPV TAENKFWRSF KVPVVAKEYS AITSIHFSAE SPYDFAVTSG ARVQIYGASS RSVKKTIA RF KDTAYSGNIR KDGKLLLAGD ATGLVQIFDL STRSILRALD AHQFPVHVTQ FCPYANTTFL SGSDDKTVKV WDLSTGAVQY DLSGHEDYVR TASWMSATRL VSGGYDGTIR LWDTRIADPE VMSFSHG EAI DVVLP MQSGS TVISAGGPSI KVWDLVAGRQ TPTK KLSNHQ KSITCLAMNS ENTRLLSGGL DGHVKIYNIS DWKVVHGMKY SGPILSMGLS PDSCNLVVG M SNGTLSQRTR RITKQTSSKT PFLGGVGS AF GVVGK KKKQIY KGENDEFYVE EARRKRLRPF DKALKSFCYS DALDMVLENG SPVLIITMLI ELLHVGG LRI ALSNRDDL SL TPLINFLRKY IRDPRFSETL CIVTSTILDI YGAALGGSVM VENAISKLRE KVEQEVAVTQ RANELVGMLQ MLSA
Source	Mammalian cell
Target Names	utp15
Protein Names	Recommended name: U3 small nucleolar RNA-associated protein 15 Short name= U3 snoRNA-associated protein 15 Alternative name(s): U3 protein 15 required for transcription
Expression Region	1-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.