



Recombinant Arabidopsis thaliana Putative flavin-containing monooxygenase YUCCA7 (YUC7)

Product Code	CSB-BP526183DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O49312
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MCNNNTSCV NISSMLQPED IFSRRCIWVN GPVIVGAGPS GLAVAADLKR QEVPFVILER ANCIASLWQN RTYDRLKLHL PKQFCQLPNL PFPEDIPEYP TKYQFIEYLE SYATHFDLRP KFNETHVQSAK YDKRFGLWRV QTVLRSELLG YCFEYICRW LVVATGENAE KVVPEFEGLE DFGGDVLHAG DYKSGERYRG KRVLVVGCGN SGMEVSLDLC NHDASPSMVV RSSVHVLPRE VLGKSTFELS VTMMKWMPVW LVDKTLVLT RLLLGNTDKY GLKRPEIGPL ELKNTAGKTP VLDIGAISMI KSGKIKIVAG IAKFGPGKVE LVDGRVLQID SVILATGYRS NVPSWLKEND LGEIGIEKNP FPKGWKGKAG LYAVGFTGGRG LSGASFDAMS VAHDIANSWK EETKQIKTV ATRHRCISH F
Source	Baculovirus
Target Names	YUC7
Protein Names	Recommended name: Putative flavin-containing monooxygenase YUCCA7 EC=1.14.13.-
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.