



# Recombinant Arabidopsis thaliana Putative flavin-containing monooxygenase YUCCA7 (YUC7)

<b>Product Code</b>	CSB-EP526183DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O49312
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MCNNNTSCV NISSMLQPED IFSRRCIWVN GPVIVGAGPS GLAVAADLKR QEVPFVILER ANCIASLWQN RTYDRLKLHL PKQFCQLPNL PFPEDIPEYP TKYQFIEYLE SYATHFDLRP KFNETVQSAK YDKRFGLWRV QTVLRSELLG YCEFEYICRW LVVATGENAE KVVPEFEGLE DFGGDVLHAG DYKSGERYRG KRVLVVGCGN SGMEVSLDLC NHDASPSMVV RSSVHVLPRE VLGKSTFELS VTMMKWMPVW LVDKTLVLVLT RLLLGNTDKY GLKRPEIGPL ELKNTAGKTP VLDIGAISMI KSGKIKIVAG IAKFGPGKVE LVDGRVLQID SVILATGYRS NVPSWLKEND LGEIGIEKNP FPKGWKKGKAG LYAVGFTGGRG LSGASFDAMS VAHDIANSWK EETKQIQIKTV ATRHRCISH F
<b>Source</b>	E.coli
<b>Target Names</b>	YUC7
<b>Protein Names</b>	Recommended name: Putative flavin-containing monooxygenase YUCCA7 EC=1.14.13.-
<b>Expression Region</b>	1-431
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.