



# Recombinant Arabidopsis thaliana Putative F-box protein At4g21240 (At4g21240)

<b>Product Code</b>	CSB-BP526190DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O49565
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDRREEEEEEE TGYGKEGTRN QSKEDESRLG KIFELIPLDM IPDILLRLPA KSAVRFRIVS KLWLSITTRP YFIRSFAPPS STRLCLMACV KARDMRLFIS LHQHDDGSYA HVDRCEIKSP KHDYYNPSSE SVNGLVCFGD FYNIVVWNPS MRQHVTLPEP KPHSTVRYFI RSCLGYDPVE DKYKVLSISG YHNGNHDPLV FTLGPQESWR VIQNSPLDIP LPTGGSRVGT CINGHVYYEA QIRFKVDDIF NFENILMSFD VRYEKFNTIK KPADPTLRNF MLNYQGKLAW FCSDYSSIRF WVLDDGDKQE WSLKNFILPF PISPEIDPIF ECRVMLTGV THTDGEFIFIN ATMCDAFCVL YYDPKSGSER RIEYEGIGDR EFCINTGILE CKRFTIDWFP NHNESLMSLV NVLKKAD
<b>Source</b>	Baculovirus
<b>Target Names</b>	At4g21240
<b>Protein Names</b>	Recommended name: Putative F-box protein At4g21240
<b>Expression Region</b>	1-417
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.