



Recombinant Arabidopsis thaliana Putative F-box protein At4g21240 (At4g21240)

Product Code	CSB-EP526190DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O49565
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDRREEEEEEE TGYGKEGTRN QSKEDESRLG KIFELIPLDM IPDILLRLPA KSAVRFRIVS KLWLSITTRP YFIRSFAPPS STRLCLMACV KARDMRLFIS LHQHDDGSYA HVDRCEIKSP KHDYYNPSSE SVNGLVCFGD FYNIVVWNPS MRQHVTLPEP KPHSTVRYFI RSCLGYDPVE DKYKVLSISG YHNGNHDPLV FTLGPQESWR VIQNSPLDIP LPTGGSRVGT CINGHVYYEA QIRFKVDDIF NFENILMSFD VRYEKFNTIK KPADPTLRNF MLNYQGKLAW FCSDYSSIRF WVLDDGDKQE WSLKNFILPF PISPEIDPIF ECRVMLTGV THTDGEFIFIN ATMCDAFCVL YYDPKSGSER RIEYEGIGDR EFCINTGILE CKRFTIDWFP NHNESLMSLV NVLKKAD
Source	E.coli
Target Names	At4g21240
Protein Names	Recommended name: Putative F-box protein At4g21240
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.