



Recombinant Feline foamy virus Protein Bel-1 (bel1)

Product Code	CSB-YP526336FEK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O56862
Product Type	Recombinant Protein
Immunogen Species	Feline foamy virus (FFV) (Feline syncytial virus)
Purity	>85% (SDS-PAGE)
Sequence	MASKYPEEGP ITEGVEEDFN SHSTSGLDLT SGNKEEPLIS LALLSMHTSK IVWIRDHFF VKILSFGGKQ KLYYICNQCH KGIPESGYIT LNTKYYLYEK GPTETGTKGL TLMRRHVQNS PCFLNSRKES GTPKTDPTRP ATSYSLCRSD YQEAGCSRPT PSNSESVCNG LGQP SERGHT SGESGGIMEE VSFISWLEGL WGEGFDYAN
Source	Yeast
Target Names	bel1
Protein Names	Recommended name: Protein Bel-1 Alternative name(s): Transactivator of spumavirus Short name= Tas Transcriptional transactivator
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.