



Recombinant *Schizosaccharomyces pombe* Uncharacterized transcriptional regulatory protein C56F2.05c (SPBC56F2.05c)

Product Code	CSB-EP526431SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60056
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDLHSILQPI STPSASVVTA PLPATIALPS PNYYYPPVAQ GHYPVNNMWS LPSNVRVISS HGAPGHQ TSA SQPSTVMPAL ETNASNAQYY PAYSVINGNN SVQVASPAYT VSHSPHSFSN PRYVAVPQKS TSPNQVCSYC EPLPNHLTKT KSCSIPPILN SDRSPLSLP TPYPVQYSTQ PVSLPQPIAA PAPPSAESSK STISDEDVAW QLIRLGALSS NSVKSSPSKS FVSISSPVQS TVKPTKASGV VKSEKVEKRS LPPQDFGNAS SSTSARRRP DHNHTSTLDA SSSNTSLAST GPMTVSSSTV ERKGKEASEV NPNSTSSVTF SDFAAAI SRS RCSRCKKSKK GCDRQRPCGR CRDAGL NSED CISDDDMPVS NARKPRGRGR GRPKTKN
Source	<i>E.coli</i>
Target Names	SPBC56F2.05c
Protein Names	Recommended name: Uncharacterized transcriptional regulatory protein C56F2.05c
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.