



# Recombinant Schizosaccharomyces pombe Homoserine O-acetyltransferase (met6)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP526432SXV   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | O60062  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MESQSPIESI VFTDSCCHPSQ QENKFVQLIS DQKIAIVPKF TLECGDILYD<br>VPVAFKTWGT LNKEGNNCLL LCHALSGSAD AGDWWGPLL<br>PGRAFDPSHF FIVCLNSLGS PYGSASPVTW NAETHSVYGP EFPLATIRDD<br>VNIHKLILQR LGVKQIAMAV GGSMGGMLVL EWAFDKEFVR SIVPISTSLR<br>HSAWCISWSE AQRQSIYSDP KFNDGYYGID DQPVSGLGAA RMSALLTYRS<br>KCSFERRFAR TVPDASRHPY PDRLPTPLTP SNAHWVVHNE GNRNRRERPC<br>RSNGSSPTSE SALNSPASSV SSLPSLGASQ TTDSSSLNQS SLLRRPANTY<br>FSAQSYLRYQ AKKFVSRFDA NCYISITKKL DTHDITRGRG SDSPKEVMKD<br>LSLPVLVLGI ESDGLFTFDE QVEIAKSFPN ATLEKIISAE GHDGFLLLEFT<br>QVNSHIQKFQ KEHLIDIMSQ TNSFERLDSQ VNDTNRESVF GEMEDITSW |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | met6  |
| <b>Protein Names</b>     | Recommended name: Homoserine O-acetyltransferase EC= 2.3.1.31<br>Alternative name(s): Homoserine O-trans-acetylase  |
| <b>Expression Region</b> | 1-489   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |