



Recombinant Schizosaccharomyces pombe Alpha-ketoglutarate-dependent dioxygenase abh1 (abh1)

Product Code	CSB-EP526433SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60066
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLAGNMEQAN VFRLEEKRYK CRADTIPDMS EVLDPNDPQS FGFEALVEIK PRVFSFQKAP GLLILKNYVS SELQMQLLKS IMFTQIQDPE NKTNLSPFYQ LPLGNDSIWR RYYNGDGESI IDGLGETKPL TVDRLVHKKL RWVTLGEQYD WTTKEYPDPS KSPGFPKDLG DFVEKVVKES TDFLHWKAEA AIVNFYSPGD TLSAHIDESE EDLTLPLISL SMGLDCIYLI GTESRSEKPS ALRLHSGDVV IMTGTSRKAF HAVPKIIPNS TPNYLLTGNK AWDGWISRKR VNFNVRQVRP SR
Source	E.coli
Target Names	abh1
Protein Names	Recommended name: Alpha-ketoglutarate-dependent dioxygenase abh1 EC=1.14.11.- Alternative name(s): Alkylated DNA repair protein alkB homolog
Expression Region	1-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.