



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C13G1.09 (SPBC13G1.09)

Product Code	CSB-EP526434SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60071
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPKAPKTKLH HAPLYKDIAE SSESGVLRQK PSKQKKSKEK NTGNGFLDAK TSKKILQLAR EQKEELDEEE NGKPSQISAF ISNGHQKDTL ENPAIESSYE ESEHARNVSD SESITSQEEE EYEELEIDDA DRDLFDRFLP TVSGEEDLT EEKTTLSLSDL IMQKINEAEA RARGEYIPSA EEEENALPPL PPKVIEVYSK VGVLLSKYRS GKIPKAFKII PTLNWDIL YLTRPDMWTP HACYEATRIF ISNLKPVQAA HFLTVILER VRDDIRENKK LNYHLYMALK KALYKPSSWF KGFLFPLVQE NCTLREAAII GSILQKVSVP VLHSAALLR LTFEFLSGAT SVFIRILLDK KYALPYKVLV SLVFYFMRWK SLERPLAVLE HQSMLVFAQR YKFDITPEQK DALLEVVRLK GHYSIGPEIR RELLSASRG EEIPVEMEY
Source	E.coli
Target Names	SPBC13G1.09
Protein Names	Recommended name: Uncharacterized protein C13G1.09
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.