



Recombinant Schizosaccharomyces pombe Phosphoglycerate kinase (pgk1)

Product Code	CSB-MP526437SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60101
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSLSTKLAIT DVDLKGKNVL IRVDFNVPLD GDRITNNARI VGALPTIKYA LEQQPKAVIL MSHLGRPNGA RVAKYSLKPV AAELSKLLGK PVKFLDDCVG PEVEKACKEA KGGEVILLEN LRFHIEEEGS AKVDGKKVKA DASAVEAFRK SLTSLGDIFV NDAFGTAHRA HSSMVGVDLP RVSGFLMKKE LDYFSKALEN PARPFLAILG GAKVADKIQL IDNLLDKVNR LIICGGMAFT FLKVLNGMKI GDSLFDDEAGS KNVESMMAKA KKNNEVEFLP VDFVTADKFD KDAKVGSA EEGIPDGWMG LDCGPKSSAK FAEVITTSKT IVWNGPAGVF EFDNFAKGTK SMLDACVKTC EAGNVVIVGG GDTATVAKKY GKEDALSHVS TGGGASLELL EGKALPGVVA LSSK
Source	Mammalian cell
Target Names	pgk1
Protein Names	Recommended name: Phosphoglycerate kinase EC= 2.7.2.3
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.