



Recombinant Putative T-box protein 32 (tbx-32)

Product Code	CSB-EP526458CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O61764
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTSQYSNYSC PLNNIQVNV I GSSRTSDDHF LEQMLTKRKK TNISQDLHYT LLGLNPRDLY QISIKLARID EKKYKFNGNF NGGDYYPHKN SEIPMEDTEE IVHLDGYQTG ARWMKTGVYF PKIGISKEQP DKERCLLLES LHAYVPVLR TTSSGISLEV TVWSQKFVTT MTSSNRSVKR REKRKLEETT GPSAKKTPEV IEFTQEAVVP EVYNAQQVYQ MQNDTMISSD SKFCSPLAFQ QDFTEEFGTI SQLDAYVNVN VQSCSSSSLS SPAALQQDST VSSDNDFDDK NSEEFDAVFE YIDQQFSNHD QNPMDQQSKS IDLNQTTWTS QSINNQS YLS TASSPAALNQ DSFASEKSSI VRDKKNDENC LDEFDRVFNE LDQQFTPTEN QLQDLSNCTT WNQEPINNSI QSNQLNNNIA PALEENVEDD NAIFEHFFIC N
Source	E.coli
Target Names	tbx-32
Protein Names	Recommended name: Putative T-box protein 32
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.