



Recombinant Arabidopsis thaliana Probable F-box protein At1g67455 (At1g67455)

Product Code	CSB-EP526494DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O64800
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MMISDLPEDM VEEILSRVSI ISLGALRWND LSKARVICKA EARQQFAGFM IKGSKVCSMR FDLHGIQNNN VEVVEPSIQ IAKFNHVEIS QVFHCDGLLL MMSKEVSNTR LVVWNPYLGK IWSIQHRSAY HSENRYTLGY DNNRKPQNLE VAYHSIKGIT YLFNGGEEKIG EGPTGCLLCF DFTRERFGPS LPLPFHVAHL YTEQLAVLLH NWGAFAMEIW VTTNIEPNAV LMEQLLKNYL LYTSDFIGSF FIDQEKHAV IFDSSLSCGS VACIIGETGN LRTVDLGVA VPYRQLPVCS YAPSLEQINQ GQELPLIGLT QIFMEQHKYW VSTFLDNLIT KYG
Source	E.coli
Target Names	At1g67455
Protein Names	Recommended name: Probable F-box protein At1g67455
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.