



# Recombinant Schizosaccharomyces pombe Uncharacterized protein C3H7.13 (SPBC3H7.13)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP526704SXV   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | O74388  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSAVITLTPL NESFQTKKLV ISPSTIYKIG RHTNKSTSPS PSNLFFNSKV<br>LSRQHAELWL DKDTLSVYIR DVKSSNGTFV NETRLSPENK PSAPCKLNSG<br>DIVDFGVDIY NEDEIVHQKV SAQVRIVVRG ATFATTPARS LDSEVMLDSI<br>MRQMFVQYQR CVELNENLES LENGVEEVSK SLSWLETGKT<br>RDNRNNHHYS RKSSPHISSL AVPSTKHLDG ERDRNLKRST SPLSSSPFVT<br>EAALNEAELA KTNLEAWKAR AFTAEARLSS KNKSWQEKKY LVLSPFFIAV<br>AGIIVYMGYW R   |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | SPBC3H7.13  |
| <b>Protein Names</b>     | Recommended name: Uncharacterized protein C3H7.13   |
| <b>Expression Region</b> | 1-301   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |