



Recombinant Schizosaccharomyces pombe tRNA-dihydrouridine (20) synthase [NAD (P)+]

Product Code	CSB-MP526734SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74731
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGLLNYSNKV CLAPMVRIGE LPMRLLALRY GANLVWGPEI VDKALLSGTP VERVVNDRIIN CIDFVKPPSN KVLFRVHPLE ANRLIFQLGS ASPELAVEAA KLVANDVAGI DLNCGCPKHF SVHAGMGAGL LKNQDRLVSI LDALVNEIGK PYKISISCKI RLLETKEDTL KLVERICDTG VRAITVHCRT TPMRNTEPAD RSYLSEIVGV CRNKDVSILV NGDVLSYNDG LDVIEKYGVD GVLIARAAER NVSCFRIEGP LSSFKVAEEF LKMALEVDNN FGNTKYCLNQ IMQGSFRKNV RQLAQTAQTY EDLKKAFEIE YKHS DASSVC PTLEKEKSLV ISFVDLPSFL ESLLASNILK QLSRLHFTKI FAVSEEEDIC KQLDDIHEKF LCHGIALSLI SADNLASAIA SAHLVICMEK DESIMNEAVC DQKITIRLPI NSNTNEAVVE CENKSMQSKH ALDIIQERIK DLEEKAQVV
Source	Mammalian cell
Target Names	dus2
Protein Names	Recommended name: tRNA-dihydrouridine(20) synthase [NAD(P)+] EC=1.3.1.91 Alternative name(s): tRNA-dihydrouridine synthase 2
Expression Region	1-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.