



# Recombinant Schizosaccharomyces pombe Set1 complex component sdc1 (sdc1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP526748SXV   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | O74861  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSNSAPARQY LNEKVTPVLL EGMKILARDR PENPLQFLGQ FLLDANANQQ<br>KQKEIVNQPE PQQETPKADA DMSTPTMAEQ VQTSFSNPAS TPLTQTSSPS<br>SNPGKNSAS   |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | sdc1  |
| <b>Protein Names</b>     | Recommended name: Set1 complex component sdc1 Short name= Set1C component sdc1 Alternative name(s): COMPASS component sdc1 Complex proteins associated with set1 protein sdc1 Lid2 complex component sdc1 Short name= Lid2C compo   |
| <b>Expression Region</b> | 1-109   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |