



# Recombinant Guanine nucleotide-binding protein alpha-11 subunit (gpa-11)

<b>Product Code</b>	CSB-EP526770CXY-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O76584
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSAADMARKN SLINRQLEKE KIDSKKMLKI LLLGGPECGK STIFKQMKII HMNGFSDL DY VNFYRLIYSN IMQSMDQLLE AAFFHFPPD DSPSIRRALN HYKSYKVRYS TSEVELNREL ADSLSKLYNA EFKSVLNRK NELKLLDSAV YFLDDIDRIS AHEYKPTMD VLRARVPTTG ITEIEFPFKQ ASLRMV DVGG QRSEQRKWIH CFDNVNGVLF IAAISGYNLY DEDEENRKDD GTPTKTNRLR YSMELFKRIA NHQCFSKKTA MILFLNKIDI FKEKIGKYPL TTCFKNYKGV NAFEPACKYV TDRFSRLVSG DIQHEKPLYT HITNATDTRN IDRVFDSCMD VIFKISMEKV GFM
<b>Source</b>	E.coli
<b>Target Names</b>	gpa-11
<b>Protein Names</b>	Recommended name: Guanine nucleotide-binding protein alpha-11 subunit
<b>Expression Region</b>	1-363
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.