



# Recombinant *Drosophila melanogaster* Protein pellino (Pli)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP526775DLU   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | O77237  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Drosophila melanogaster</i> (Fruit fly)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MVKRTDGTES PILAEDGGDG HDKPRLRYGE LVILGYNGYL PQGDRGRRRS<br>KFVLHKRTEA SGVKRSKHYI VQSPQTSKAI LDANQHSISY TLSRNQAVIV<br>EYKEDTETDM FQVGRSSESP IDFVVM DTL P GDKKDAKVMQ STISRFACRI<br>LVNRCEPAKA RIFAAGFDSS RNIFLGEKAT KWQDNVEIDG LTTNGVLMH<br>PKGSFCGGNA KCGLWRECSV GGDVFSLRES RSAQQKGQPI YDECNILQDG<br>TLIDLCGATL LWRSAEGLQH SPTKHDLEKL IDAINAGRPQ CPVGLNLT LVI<br>PRKVNIGDQV NQPYVYL NCG HVQGHHDWGQ DENTGARRCP<br>MCLELGPVVT LCMGLEPAFY VDVGAPTYAF NPCGHMATEK TVKYWANVEI<br>PHGTNGFQAV CPF CATPLDG ATGYIKLIFQ DNLD |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | Pli   |
| <b>Protein Names</b>     | Recommended name: Protein pellino   |
| <b>Expression Region</b> | 1-424   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |