



Recombinant Arabidopsis thaliana Probable protein phosphatase 2C 5 (At1g09160)

Product Code	CSB-YP526808DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O80492
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSVSKASRTQ HSLVPLATLI GRELRSEKVE KPFVKYQGAA LAKKGEDYFL IKTDCERVPG DPSSAFSVFG IFDGHNGNSA AIYTKHELLE NVVSAIPQGA SRDEWLQALP RALVAGFVKT DIEFQQKGET SGTTVTFVII DGWTITVASV GDSRCILDTQ GGVVSLLTVD HRLEENVEER ERITASGGEV GRLNVFGGNE VGPLRCWPGG LCLRSIGDT DVGEFIVPIP HVKQVKLPDA GGRLIISDG IWDILSSDVA AKACRGLSAD LAAKLVVKEA LRTKGLKDDT TCVVVDIVPS GHLSLAPAM KQNPFTSFL SRKNHMDTNN KNGNKL SAVG VVEELFEEGS AVLADRLGKD LLSNTETGLL KCAVCQIDES PSEDLSSNGG SIISSASKRW EGPFLCTICK KKKDAMEGKR PSKGSVTT
Source	Yeast
Target Names	At1g09160
Protein Names	Recommended name: Probable protein phosphatase 2C 5 Short name= AtPP2C05 EC= 3.1.3.16
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.