



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At2g44030 (At2g44030)

Product Code	CSB-BP526811DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O80573
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTTEEDMRSS SSESSPKSFL SLPYDVVFNC LSRVSRTHDP ILSLVSKSFR SLLALPDLEA ERFRILKNET CLYVCLNLNN NNNPNPSWFI LSQTPKHKLI PLPSLPYPDP HPNCSTVVST GSEIYLLGGF VAKEKRSRRA YVLDCKSHQW RRLPKMRIAR KEAAANVIDG KINVYGGCSS EYHNSVNWGE IYDPMTQTWE PFPEGALNKE GVIPCALIKD GIAFPDCGLL ISGKVYDTTT MDTLDYNTPM DKLDLCPNVC MLKIDNQDFQ ASVSDGKCLKL VRCRGAMAWH WTVGGLEELS CNYLLSVASP GGGRRVTVWW KTYTKECKTE IWCALILLER ELVWGVIEWS ENVFTLPGSE SDSDSNSFLL HYELVTLARH
Source	Baculovirus
Target Names	At2g44030
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At2g44030
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.