



# Recombinant Arabidopsis thaliana Heat stress transcription factor A-1b (HSFA1B)

<b>Product Code</b>	CSB-YP526843DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O81821
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MESVPESVPS PNSNTPSIPP PVNSVPPFLS KTYDMVDDPL TNEVVSWSWG NNSFVWWSAP EFSKVLLPKY FKHNFFSSFV RQLNTYGFRK VDPDRWEFAN EGFLRGRKQL LKSIVRRKPS HVQQNQQQTQ VQSSSVGACV EVGKFGIEEE VERLKRDKNV LMQELVRLRQ QQQATENQLQ NVGQKVQVME QRQQQMMSFL AKAVQSPGFL NQLVQQNNND GNRQIPGSNK KRRLPVDEQE NRGDNVANGL NRQIVRYQPS INEAAQNMLR QFLNTSTSPR YESVSNPDS FLLGDVPSST SVDNGNPSSR VSGVTLAEFS PNTVQSATNQ VPEASLAHHP QAGLVQPNIG QSPAQGAAPA DSWSPEFDLV GCETDSGECF DPIMAVLDES EGDAISPEGE GKMNELLEGV PKLPGIQDPF WEQFFSVELP AIADTDDILS GSVENNDLVL EQEPNEWTRN EQQM KYLTEQ MGLLSSEAQR K
<b>Source</b>	Yeast
<b>Target Names</b>	HSFA1B
<b>Protein Names</b>	Recommended name: Heat stress transcription factor A-1b Short name= AtHsfA1b Alternative name(s): AtHsf-18 Heat shock factor protein 3 Short name= HSF 3 Heat shock transcription factor 3 Short name= HSTF 3
<b>Expression Region</b>	1-481
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.