



Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 20 (At2g19810)

Product Code	CSB-MP526857DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82199
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MMIGESHGRGF NPTVHIPWP LSEDLTVSDI YGSPDGGSSM MEALAEALQRY LPSNEPDPDS DPDLSPDPS IDAYTCDHFR MYEFKVRRCR RGRSHDWTEC PYAHPGEKAR RRDPRKFHYS GTACPEFRKG CCKRGDACEF SHGVFECWLH PARYRTQPCK DGGNCRRRVC FFAHSPDQIR VLPNQSPDRV DSFDVLSPTI RRAFQFSISP SSNSPPVSPR GDSDSSCSLL SRSLGSLGN DVVASLRNLQ LNKVKSSLSS SYNQIGGYG SGFGSPRGSV LGPGRSLPT TPTRPGFMNI WENGLEEPA MERVESGREL RAQLFEKLSK ENCMGRIEPD PDQGAGDTPD VGWVSDLVM
Source	Mammalian cell
Target Names	At2g19810
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 20 Short name= AtC3H20
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.