



Recombinant Arabidopsis thaliana Probable transcriptional regulatory protein At2g25830 (At2g25830)

Product Code	CSB-MP526864DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82314
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASHCSMRVI LLRFSNGVSS RSILNSTNHR LLSLTMTNTL SSLSSISPHT TTSHF TAASQ QSDDQNCFRK LQLRKISIST PLCMGRRSSK IAGRKGAQDS KKAKLYCRIG KEVVS AVKKG GPNPVSNTTL ATILDKAKEL DVPKDIVERN IKRASEKGQE AFIEKIYEVY GYGGVSMVVE VLTDKINRSV AAIRSVVKDY GGMADSGSV MFKFKRVRV V NIKVTEADKD QLLI ALDAG AEDVIEPPTY EDDTDEDREE RYYKIVTSNE NYSTILSKLR DEGVNFEPDN GSELLPLTTV EVDDEAMELN KELMQKLEL DDVDAVYIDQ K
Source	Mammalian cell
Target Names	At2g25830
Protein Names	Recommended name: Probable transcriptional regulatory protein At2g25830
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.