



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At2g29810 (At2g29810)

Product Code	CSB-BP526867DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82375
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAVISETSDE GSNGGAANKK PQEEEEENIPP IPKELPEELI VIIVALVRRY HYPKLSLISK AYRDLISSPE LFQTRSRLGF TEPVLYTSIG FPPFDLPSWY ILHRISLQFK QITSLPSMLP GSAVVTIDYK MYVLGGFIGL NQPVSTMIVI DCRFHTYREL PSMQRDRGGA AAGVIDGKIY VIGGCKKRYN DWVEVFDVEN ESWETVPGPY PNVASESVEF SQYAVMEQKI YILDPLFCLT YEPRKRRFEA WGRRSQRRY ESCCVNDML CTVDSTCSHG YPIVVYDPKG MVWRPVKGVQ SSDLPNLVYY ESRMANFGGK LVILGGNQSR DRKDSYLEKD IWCIEIALEK REDGDIWGHV ESRSVVLSSS KWPLIALSRS VTV
Source	Baculovirus
Target Names	At2g29810
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At2g29810
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.