



Recombinant Arabidopsis thaliana Rac-like GTP-binding protein ARAC10 (ARAC10)

Product Code	CSB-EP526871DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82481
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASSASKFIK CVTVGDGAVG KTCMLICYTS NKFPDYIPT VFDNFSANVV VEGTTVNLGL WDTAGQEDYN RLRPLSYRGA DVFVLSFSLV SRASYENVFK KWIPELQHFA PGVPLVLVGT KLDLREDKHY LADHPGLSPV TTAQGEELRK LIGATYYIEC SSKTQQNVKA VFDSAIKEVI KPLVKQKEKT KKKKKQKSNH GCLSNVLCGR IVTRH
Source	E.coli
Target Names	ARAC10
Protein Names	Recommended name: Rac-like GTP-binding protein ARAC10 Alternative name(s): GTPase protein ROP11
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.