



Recombinant Arabidopsis thaliana Putative F-box protein At4g10740 (At4g10740)

Product Code	CSB-EP526872DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82488
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRTTMSNLPK ELVEDIVSRV PLHCLRAMRL TCKNWNALLE SQSFKKMHIR KEEEAARELG ENRMIVMMDY NVYLMGIIVN ENLSIKSLGK LTCLHGSKQV KISQVFHCEG LLLCILKDDD TKIVVLNPYL GQTRWIQTRK YHRTSEWEGR DVYNYALRYE TNSGNRSPKI LRFIDDFHHH PENPALRYET YDFDSDLWTT LDVSPHWRII SRCGGVSLNG NTYWGAVERN ASAYIGHMIC FDFTRGGFGP LIPLPFKARG LVFAPLSSDR AEKIAALFRT SKKVVIWITT KIDATNVTWS NFFTIIYISNL ERNLSCKSFF TIYISNLERN LSCKSFFIDM EKKVAMVFDK EKGKKVAHNI INIIGEAGCV RKLVLKESGD KSCWPLVCSY VPSTVQIKQH NRGKRQKKVI RKGIDIDMRH DKRIKVD
Source	E.coli
Target Names	At4g10740
Protein Names	Recommended name: Putative F-box protein At4g10740
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.