



Recombinant Schizosaccharomyces pombe Probable U3 small nucleolar RNA-associated protein 16 (utp16)

Product Code	CSB-MP527172SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94259
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MQSSRKIMIL KRNIKGLKSS AIIATYLTLN KLFIDKFCKT QESTYPIHFV HCHFGGIKKF HSHCQLSSSY TNSGQKSFKS FNSLRFIVLK SKTCLLRYKL LIKTFCCITY SKIENLSNTR MARLTAPSKG PSFSNLNKPKN DLLKAEIEDN EASHYSQNQS SDSDSDEAPE EVSLKSSSDE VKKRLQRLKG NELQLRKAER EKRRLLQNERN RERNAQRVSA SKLDLAILEA AEVEEINSTT EIEEKVDEAS LDLYTPKGHK IVFPNNEIPK KRTRSFKRGD FKVSVLKETN DPLLAPNYEK KLSQRKRDWI NRQSVPRASK RRRPLTY
Source	Mammalian cell
Target Names	utp16
Protein Names	Recommended name: Probable U3 small nucleolar RNA-associated protein 16 Short name= U3 snoRNA-associated protein 16
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.