



Recombinant *Schizosaccharomyces pombe* Cystathionine beta-lyase (str3)

Product Code	CSB-MP527189SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94350
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MPSDCKYSVD TELVHVEGNE DQYHASSVPI YQSATFKQPC LEHMGKFDYT RSGNPTRSVL QVHLAKLMKA KHAFVTSNGM SALDMILRCC KSNSHVVAGH DLYGGSDRLL SFNQRQYGFK VDNVDTSDLA AFEAALRPDT NLVLIESPTN PRISICDIRA IVKITRSKAK DALLVMDNTM LSPVLCNPLD FGYDIVYESA TKYLSGHHDL MGGVIATKSD EIAKSVFFNI NAMGAAMAPF ECFLLLRGIK TMGLRVERAQ QNAIEIAKFL KSKGLQVFP GLDPDAKSTA IFYSFARGPG AVMSVFTGDV EVSKTIVNTT KLFEISVSFG AVNSLISMPA YMSHASIKKE VRDARGLSED LIRICVGIEN VDDLKADLEN ALAQANFKQN
Source	Mammalian cell
Target Names	str3
Protein Names	Recommended name: Cystathionine beta-lyase Short name= CBL EC= 4.4.1.8 Alternative name(s): Beta-cystathionase Cysteine lyase
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.