



Recombinant *Schizosaccharomyces pombe* Uncharacterized RNA-binding protein C660.15 (SPBC660.15)

Product Code	CSB-YP527205SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94432
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGGSDFEDDE LFKDLYGEEN EKKVESASGN QETSNVTPTK ENEGYEELEK SGEAGAERTK ENPFREEPGA DFDRSGSDQQ YSAEAEANQE DDLNETSQQA PDPSSYDIRG QALSSATWDN AEDGENSKND NYNENQSALT GSGAMESNED NAEETSPFNR EDGKMFIGGL NWETTDDSLR DYFEQFGEVL DCTVMRDSTT GRSRFGFGLT FKNPKCVNEV MSKEHHLDGK IIDPKRAIPR EEQEKTAJMF VGGVPGDCTE EEFNRFFNQF GRVLDATLMM DKDTGRPRGF GFVTYENESA VEATMSQPYI TIHGKPVVEVK RATPKASLRD SHDRHQHGYH GNANPYAQN MNMYGGMTPA MMAQYYRQMQ QYMEAMRNMP AAAGAVPYPQ PVMPAGMADW QQQQQQGAAY FDPSKMNQGT GDGVFPSPSM PSGSSRGGYH GRNPGGPNRQ RYRGRRDGRG GNTGGGHSFH PYRR
Source	Yeast
Target Names	SPBC660.15
Protein Names	Recommended name: Uncharacterized RNA-binding protein C660.15
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.