



Recombinant Schizosaccharomyces pombe Uncharacterized PH domain-containing protein C637.13c (SPAC637.13c)

Product Code	CSB-YP527207SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94447
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MELSGKSEFP TRTEIPNQAS NGDAVPNTDA YMGLCNSLTY RFEGWRHLVE NLIAFFKQLE STSKNNAKEY MKLSKVIVHP YKDANVFEEH GIQDIFAALR DQTAAIASDS EQISIQLPGT IINILELLRE DLREHCKKIA AEGTKGVKVV EKQRAETQKY LSLLD RALLP WRASSPPTVN IKNDPFIVDR QVLNCLARQV QEENNHSVAV AQLQEFSFRF EQNLIKIKD TVKQFEGMMN QTHVKAINHL QEVVRVSEAQ TLAGEWGTYA KREPEFIHGS VPPRSVDAIK YPGKNDQPTV PIMAGYLIRK TSFLKKKQRG FYAFTHSGYL YEFKSSDSLQ DPEPEFALYI PDCLIGRPSE KKP KFKITGK DATKKIFGAR NDYAFRASNN TELMRWWEAL NTHIANVNYI QPLSNANGPS AVSDSDDDDD DPNDFRPAVE RQSSTMNTRM SQPSSAVNTN RSYGSEQIPS YADSQGN SGA NNNFIYNASA TGHNAWNV
Source	Yeast
Target Names	slm1
Protein Names	Recommended name: Uncharacterized PH domain-containing protein C637.13c
Expression Region	1-498
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.