



# Recombinant *Drosophila simulans* Protein anon-37Cs (anon-37Cs), partial

<b>Product Code</b>	CSB-BP527271DMJ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O96566
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila simulans</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FKLASRRRLY NARVLQADNI GDKQRNPDVD AARQNTQIVV VGAGLPGLSA AQHLLYNGFR RTVILEATDR YGGRINTQRF GDTYCELGAK WVKIDGSQDS MYELLRNTEG LDKQIKQPDR ATYLQDGS HI NPAMVELIDT LFRQLCRGFK VSERVKTGGD LHSLDNVMNY FRTESDRIIG TSFQQPKDQL AAREIFQSLF KEFGSILGCC LEYVNIEHIT KCPVQQELRP LYVPTGLDNV VDDLIQNMDK AQLQTGKPVG QIQWTPAPMK SVGCLDGS LY NADHIICTLP LGVLKSFAGV LFRPTLPLDK MLAIRNLGFG NPLKIYLSYK KPIGRWLKGS LRPLGTLNLP SAEQQPERNW TQQVVEISQV PSSQHVLE VH VGGGYEEIE KLPDDELLEQ ITGLLRRCVS NNLVPYPQEL LRSNWSTSAC YLGGRPYFST INSARDVQRL AAPLGEKSPG L LFAGDATSL NGFGTIDAAR SSGIREAQRI IDFYLKRAHF G
<b>Source</b>	Baculovirus
<b>Target Names</b>	anon-37Cs
<b>Protein Names</b>	Recommended name: Protein anon-37Cs
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.