



Recombinant Schizosaccharomyces pombe

Uncharacterized protein C16C6.04

(SPBC16C6.04)

Product Code	CSB-YP527311SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42928
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDSFHPTPGK PTTATSNSSL NFFVRNREKS IRCSSLHQTF TAPSPEASPS IGLRYVNYSM SDETDESMFP LDELTDMED DDTEYISDST TDLPSAAMAR GTRQLSYNEN LDQRSFSKTP NKHIEVKNLK DLCSPSHSGR ISKSLCPVLR RKSLLPKPKM FQRVASALYE ESSPLELEIR SESEFSKMFF EPKSSSFVS RAASPAMVGP GVPFYSSITG ANGNVLAPSG SLKAVENS DV IESSAEDSSN TDNPSTKPSN EMPISPLLSS SVFFKDTEMS TPNSNHSRSR TPSSKKRTRW GEEIIDLSKR RAVSPSIYYD LDKKCSPIHS VMVSPKIKD THEVLMNLKL
Source	Yeast
Target Names	SPBC16C6.04
Protein Names	Recommended name: Uncharacterized protein C16C6.04
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.