



Recombinant Schizosaccharomyces pombe Protein shq1 (shq1)

Product Code	CSB-MP527321SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43076
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MITPKFELRQ DEEFVYLDVH TPYIKAKNVE IDAFNNEVNF CIAPYFLKLE LPGNVLDDDR ANASYDISSG FVHIKFPKET PGEVFPDLDL LTTLLVNQNK EPVKKPLIEE VEDTPSPSTS QNASKDEEKE VKIDWQNIHF DDPDFDWGLL QRLDDEKMQY TSTSQYGFNN QYSGLLK YNA LVGNEINQIP EPERTPPSTR SEIRVQLEDE KFDAEHYMA D FYDREMIDEI LHYIPDYIGE FLLVQQGKSD RLLEFTSEEK EALAKLPKRS YLIDNPKNVY LCLVPLIFAY AYDNRTTLGD PTIESCW NIG TLSSTISCLD TNFLAIKDIF IACIRRS LAY PLYRNWDLAL TCWKDTHAIL SLGKRWILRV LLRIRNLFEH HDVYYIYNAI VLNDYAYWIQ TANDNVLAAL AYQVSECTIT KEETNWPLEE IEANVQETII ASDEESESIE K
Source	Mammalian cell
Target Names	shq1
Protein Names	Recommended name: Protein shq1
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.