



# Recombinant *Schizosaccharomyces pombe* Zinc-type alcohol dehydrogenase-like protein C16A3.02c (SPBC16A3.02c)

<b>Product Code</b>	CSB-MP527335SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O42909
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSETKSEFG TKSALYLNRT GKPKDVLYLE KGLHIPNPAE LGPYDVLVEV VATSINPLDY KLMNTYQMIA KALFKLPNIP GYDFAGRVLA VGSEVKEFSA TQRVWGCQSF PRAGRQGGSC ATHIVTGDKD VWHLDPDGVSF NEGAGFGIAG LTAWEVLRQ MKVKPGTKLV IEGASGGVGT FAVALAKALE CEVTTISSTE NLDLCKSLGA THTLDYKKDN LVERLADLGP YDFVFDCVND NVLYRASSKF VKPDGAFFGI GGDITLSYVG SRLSRTLPRP VLGSSSHSYY NILLHVDQEM LRDFVDFVMK HNIKTVIDSV YDFEDTVEAF NRLMTHRCKG KVIKTD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SPBC16A3.02c
<b>Protein Names</b>	Recommended name: Zinc-type alcohol dehydrogenase-like protein C16A3.02c EC= 1.-.-.
<b>Expression Region</b>	1-347
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.