



Recombinant *Paeonia suffruticosa* Maturase K (matK)

Product Code	CSB-MP527382EQO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O46993
Product Type	Recombinant Protein
Immunogen Species	<i>Paeonia suffruticosa</i> (Tree peony) (<i>Paeonia moutan</i>)
Purity	>85% (SDS-PAGE)
Sequence	MEKSQGYLEL DKSWRHDFLY PLIFQEYIYA LAHEQGLNRS ILLENTDHDN KYSSLIVKRL ITRMQQQNHF LIFDNDSNQN PFWKHNNNLY SQMISEGFVI IVEIPFSPRF VDSLEEKKKI LKFNNLRSIH SIFPFLEDQI LHLNFVSNIL IPYPIHLEIV VQSLRYRVKD ASSLHLLRFF LFTLNKSISS FSKRNPRLFL FLYNSHVY EY ESTFLFLRNK TSHLRSTSSG AFLERIFFYG KIKHLIEVFA NDFQAILWLF KDPFMHYVRY QGKSILASKR TSLRMNKWKY YLINFWQCQF YVWSQPGRVS INQLSNHSLD FLGYLSSVRL NPLAVRSQML ENSFLT DNAI KKFDIIVLLI PLIGSLAKAK FCNVLGHPLS KPARADSSDS DIERFVRIC RNLSHYHSGS SKKKS LYRIK YILRLSCART LARKHKTA VR SFLKRLGSEL LEEFLT EDGQ VISLIFPRTS STSWRLYRGG IWYLDITCIN DLANHE
Source	Mammalian cell
Target Names	matK
Protein Names	Recommended name: Maturase K Alternative name(s): Intron maturase
Expression Region	1-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.