



Recombinant Arabidopsis thaliana Putative F-box protein At2g33190 (At2g33190)

Product Code	CSB-EP527406DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O49316
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSRDDNGWSK LYPDLLRSIF ESLSCLDFHR AGTVCSNWYA VSRSCPLYPW RIVFRGKNSV LFDPIQDKIY TKNLLGIDLS KIHCLASYGN WILIVDPRLD FYLLNVFTRE TINLPSLESS LRGDRPFRFI RNDDESGYFL ELYGTKFLLT WNDFRFEKTA ILWVNGRNGD YVVAWAIKQF YIFSYYKIGN DDDDKRWSIT CTQCEDMAYK DNKLYVYTFD HYINILDFSG NSPKPEMEEN PYLSHPFSFV DAIYKLRLAV TREGKVLIVL SLVGLDKKIC IYKLNKIGD WEIVESLGDE MLIFGHGVTI RAPDKDISGG GLKSDSICFL YDDYLSDHYL TMKRQPICGV FDLATSTITW HTRLEDLSSK ICWFVPGCA
Source	E.coli
Target Names	At2g33190
Protein Names	Recommended name: Putative F-box protein At2g33190
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.