



# Recombinant Leptospira interrogans serogroup Icterohaemorrhagiae serovar Lai Flagellar filament 35 kDa core protein (flaB)

<b>Product Code</b>	CSB-BP527481LOS
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O51941
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Leptospira interrogans serogroup Icterohaemorrhagiae serovar Lai (strain 56601)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIINHNL SAV NAHRSLKFNE LAVDKTMKAL SSGMRINSAA DDASGLAVSE KLRTQINGLR QAERNTE DGM SFIQTAEGFL EQTSNIIQRI RVLAIQTSNG IYSNEDRQLV QVEVSALVDE VDRIASQAEF NKFKLFEGQF ARGSRVASMW FHMGPQNQR ERFYIGTMTS KALKLVKADG RPIAISSPGE ANDVIGLADA ALTKIMKQRA DMGAYYNRLE YTAKGLMGAY ENMQASESRI RDADMAEEVV SLTTKQILVQ SGTAMLAQAN MKPNSVLKLL QQI
<b>Source</b>	Baculovirus
<b>Target Names</b>	flaB
<b>Protein Names</b>	Recommended name: Flagellar filament 35 kDa core protein Alternative name(s): 35 kDa antigen Flagellin class B
<b>Expression Region</b>	1-283
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.