



Recombinant Human Exocyst complex component 3 (EXOC3)

Product Code	CSB-YP527663HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60645
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MQCEDSTSFF TMKETDREAV ATAVQRVAGM LQRPDQLDKV EQYRRREARK KASVEARLKA AIQSQLDGVR TGLSQLHNAL NDVKDIQQSL ADVSKDWRQS INTIESLKDV KDAVVQHSQL AAAVENLKNI FSVPEIVRET QDLIEQGALL QAHRKLMdle CSRdGLMYEQ YRMDSGNTRD MTLIHGYFGS TQGLSDELAK QLWMLQRSL VTVRRDPTLL VSVVRIIERE EKIDRRILDR KKQTGFVPPG RPKNWKEKMF TILERTVTTR IEGTQADTRE SDKMWLVRHL EIIRKYVLDD LIVAKNLMVQ CFPPHYEIFK NLLNMYHQAL STRMQDLASE DLEANEIVSL LTWVLNTYTS TEMMRNELA PEVDVGTLEP LLSPHVVSEL LDTYMSTLTS NIIAWLRKAL ETDKKDWWKE TEPEADQDGY YQTTLPAIVF QMFEQNLQVA AQISEDlKTK VLVLCCLQQMN SFLSRYKDEA QLYKEEHLRN RQHPhCYVQY MIAIINNCQT FKESIVSLKR KYLKNEVEEG VSPSQPSMDG ILDAIAKEGC SGLLEEVFLD LEQHLNELMT KKWLLGSNAV DIICVTVEDY FNDFAKIKKP YKKRMTAEAH RRVVVEYLRA VMQKRISFRS PEERKEGAEK MVREAQLRF LFRKLASGFG EDVDGYCDTI VAVAEVIKLT DPSLLYLEVS TLVSKYPDIR DDHIGALLAV RGDASRDMKQ TIMETLEQGP AQASPSYVPL FKDIVVPSLN VAKLLK </p>
Source	Yeast
Target Names	EXOC3
Protein Names	Recommended name: Exocyst complex component 3 Alternative name(s): Exocyst complex component Sec6
Expression Region	1-756
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.